

Implant Compatibility & OI Part #s

Astra Tech® Osseo Speed

Platforms	3.5/4.0	4.5/5.0
Titanium-Base Set (engaging)	GH 1.15, H 10 3721	GH 1.15, H 10 3722
Titanium-Base Set (non-engaging)	GH 1.15, H 10 3750	GH 0.8, H 10 3751
Screw	3769	3770
Analog	3614	3615
Scan Body	3664	3665

Key:
Titanium (Ti) Base Sets include a screw.
GH= Gingival Height (mm)
H= Chimney or Post Height (mm)

Libraries available at:
openimplants.com/downloads

Biomet® 3i Certain


Platform	3.4	4.1	5.0
Titanium-Base Set (engaging)	GH 0.8, H 10 3707	GH 0.5, H 10 3708	GH 0.5, H 10 3709
Titanium-Base Set (non-engaging)	GH 1.5, H 8 5813	GH 1.5, H 8 5814	
Titanium-Base Set (non-engaging)	GH 0.8, H 10 3736	GH 0.5, H 10 3737	GH 0.5, H 10 3738
Screw	3763	3763	3763
Analog	3598	3599	3600
Scan Body	3652	3653	3653

Nobel Active®

Platform	3.5 NP	4.3/5.0 RP
Titanium-Base Set (engaging)	GH 0.8, H 10 3702	GH 0.5, H 10 3703
Titanium-Base Set (engaging)	GH 1.5, H 8 5815	GH 1.5, H 8 5817
Titanium-Base Set (engaging)	GH 3.0, H 8 5816	GH 3.0, H 8 5818
Titanium-Base Set (non-engaging)	GH 0.8, H 10 3731	GH 0.5, H 10 3732
Screw	3758	3759
Analog	3593	3594
Scan Body	3647	3648

Nobel Biocare® Multi-Unit Compatible

Titanium Insert	Ti Temporary Cylinder	Scan Body	Analog	Prosthetic Screw
-----------------	-----------------------	-----------	--------	------------------



Nobel® Replace

Platform	3.5 NP	4.3 RP	5.0 WP	6.0
Titanium-Base Set (engaging)	GH 0.8, H 10 3698	GH 0.5, H 10 3699	GH 0.5, H 10 3700	GH 0.5, H 10 3701
Titanium-Base Set (engaging)	GH 2.0, H 8 5819	GH 2.0, H 8 5820		
Titanium-Base Set (non-engaging)	GH 0.8, H 10 3727	GH 0.5, H 10 3728	GH 0.5, H 10 3729	GH 0.5, H 10 3730
Screw	3756	3757	3757	3757
Analog	3589	3590	3591	3592
Scan Body	3643	3644	3645	3646

Straumann® Bone Level

Platform	3.3 NC	4.1/4.8 RC
Titanium-Base Set (engaging)	GH 0.9, H 10 3713	GH 1.0, H 10 3714
Titanium-Base Set (engaging)	GH 2.0, H 8 5821	GH 3.0, H 8 5823
Titanium-Base Set (engaging)	GH 3.0, H 8 5822	
Titanium-Base Set (non-engaging)	GH 0.8, H 10 3742	GH 0.8, H 10 3743
Screw	3765	3765
Analog	3604	3605
Scan Body	3656	3657

Zimmer® Tapered Screw Vent

Platform	3.5	4.5	5.7
Titanium-Base Set (engaging)	GH 0.85, H 10 3718	GH 0.25, H 10 3719	GH 0.25, H 10 3720
Titanium-Base Set (engaging)	GH 1.5, H 8 5824	GH 1.5, H 8 5826	
Titanium-Base Set (engaging)	GH 3.0, H 8 5825		
Titanium-Base Set (non-engaging)	GH 0.8, H 10 3747	GH 0.25, H 10 3748	GH 0.25, H 10 3749
Screw	3768	3768	3768
Analog	3611	3612	3613
Scan Body	3661	3662	3663



Learn more at openimplants.com

Order online or call toll free: 1-888-800-8369